

# Software

## Update D-Lab 3.0, 3.01 or 3.10

Before you update you D-Lab 3.0, 3.01 or 3.10 to the latest D-Lab version, it is recommended to [backup your data](#) and additionally the old D-Lab version.

Therefore make a copy of the folder C:\Programm\D-Lab and rename this copy to “D-Lab 3.xx”. You will then have the option to downgrade to your old D-Lab version.

After this backup download the [new D-Lab 3.xx Setup](#) and install it.

Unique solution ID: #1186

Author: Ergoneers Gmbh

Last update: 2017-04-27 14:13