

# Getting Started

## How to install D-Lab?

Before you download D-Lab please check the [hardware requirements](#).

Then download the latest D-Lab setup [here](#) and follow the installation wizard. To run the software a valid licence dongle is required.

In case you already installed a previous version make sure you have a [backup](#) of all your data. Afterwards follow the same procedure as for the first installation.

Updating a D-Lab 3.0, 3.01 or 3.10. please this additional information [here](#).

Unique solution ID: #1007

Author: Ergoneers GmbH

Last update: 2017-04-27 14:14